


# Design meets function Knee joint series KINEGEN



# KINEGEN series - sophisticated prosthetic knee joints with high functionality

Prosthetic knee joints are intended to replace the function of a lost limb and assist in the restoration of an amputees ability to walk. In order to achieve this, the device needs to meet high technical and quality standards. This ensures that the patient can lead an independent and carefree life. Which knee joint is used depends primarily on the mobility and weight class of the patient. For this reason, close attention must be paid to the precise identification of these criteria to ensure best possible care for the patient.

Prosthetic knee joints not only have to replace the biological knee joint, but in many cases replace the functionality of a large part of the entire missing limb. The prosthetic knee joint must meet numerous requirements - it should provide the patient with a high degree of safety, correspond to the dynamics of the prosthesis wearer and be adaptable to the patient's gait pattern.

The KINEGEN series therefore offers an excellent selection of technically sophisticated prosthetic knee joints for every need.

The high functionality combined with an aesthetic design enables the prosthesis wearer to enjoy a better quality of life, both in everyday life and in leisure time

## The KINEGEN series at a glance



#### 3A860 KINEGEN.smart

Stability is especially important for geriatric walkers. Locked joints are often used, which release the flexion of the lower leg (e.g. for sitting) only after unlocking.

These demands are taken into account in the KINEGEN.smart, which has a load-dependent brake and a lock that can be deactivated.



### 3A1000 KINEGEN.air

Swing phase-controlled joints for medium walking speeds such as the KINEGEN.air are equipped with a load-dependent brake, knee-extending spring element and pneumatics. This allows for continuous adjustment of the movement resistances for flexion and extension as well as the end position damping. This comfortably dampens the extension stop just before full extension of the joint.



## 3A900 | 3A900/1 KINEGEN.wave

The KINEGEN.wave is a polycentric, water-proof knee joint. The pivot point is offset to the rear of the knee joint which increases the stance phase stability. The safety is futher enhanced through a lock which can be deactivated for higher mobility patients. For the first time, it is possible to bathe, shower or swim with a knee joint from the KINEGEN series.



### 3A1800 KINEGEN.air-active

The KINEGEN.air-active is a polycentric knee joint with pneumatic swing phase control. The polycentrics allow for a particularly large flexion angle and offer more trip safety due to the leg shortening in the swing phase.



## 3A2500 | 3A2500/4 | 3A2500/3 KINEGEN.stream

The polycentric knee joint KINEGEN.stream is designed for highly active patients. Due to the powerful hydraulics, an optimal individual adjustment to the prosthesis wearer is possible. Various connections are possible proximally (pyramid and M36 thread).

## 3A860 KINEGEN.smart Brake Knee Joint with Lock

With its colour and shape concept, the KINEGEN.smart 3A860 fits harmoniously into the knee joint series. Despite the particularly low mounting height and the low own weight, it achieves a weight approval of up to 125 kg and enables the technician to cover a particularly wide range of patients.

If the integrated lock is released for maximum stability, the KINEGEN.smart can be used as a brake knee joint. Here, both the braking force and the extension assist can be individually adjusted to the user.



















# 3A1000 KINEGEN air Brake Knee Joint with Pneumatic Swing Phase Control

A fundamental benefit of the KINEGEN.air is the precisely adjustable shock absorbing action which damps the extension impact comfortably shortly before complete extension. Therefore the patient will not experience exposure to heavy or unnecessary forces on the residual limb and the gait pattern will be considerably smoother. Depending on the desired adjustment it is possible to achieve a completely impactfree feeling in the residual limb during knee-joint extension. High stance stability through load-dependent brake and posterior positioned joint axis.















adjustable hydraulic end position damping for optimal wearing comfort

Simple regulation of the pneumatics - can be adjusted independently of each other via flexion and extension damping



# 3A900 KINEGEN.wave Polycentric Waterproof Knee Joint with Lock

A major plus point of the KINEGEN.wave 3A900 is its water resistance: Prosthesis wearers can move carefree and freely in fresh, salt or even chlorinated water with the polycentric knee joint. This increases flexibility and comfort for the user. Geriatric walkers in particular also require a maximum safety standard, which the KINEGEN.wave optimally fulfils with a high stance phase safety and a lock that can be deactivated. The robust joint is the ideal interim and initial fitting and is perfect for use in rehabilitation therapies in the swimming pool.










The KINEGEN.wave is available in two versions. The following proximal connections are available:



Adjustment core (3A900)



M36 thread (3A900/1)



# 3A1800 KINEGEN.air-active Polycentric Knee Joint with Pneumatic Swing Phase Control and Integrated Extension Assist

The polycentrics allow for a particularly large flexion angle and offer more trip safety due to the leg shortening in the swing phase. The pneumatic swing phase control in combination with the dynamics of the multi-joint chain ensures a harmonious and natural gait pattern for the prosthesis user. The pneumatic extension and flexion resistance can be continuously adjusted. The KINEGEN.air-active is low-maintenance and intuitive to operate.



















# 3A2500 | 3A2500/3 | 3A2500/4 KINEGEN.stream Polycentric Knee Joint with Hydraulic Swing Phase Control

The KINEGEN.stream polycentric knee joint combines all the advantages for highly active users. The axle arrangement allows for maximum dynamics, while a powerful hydraulic unit with precise adjustment options enables individually optimized adaptation of the joint to the prosthesis user. The infinitely adjustable end-position cushioning helps to increase wearing comfort.













3A2500/3 (as 3-Arm Socket Adapter is approved up to max. 125 kg body weight) The KINEGEN.stream is equipped with a free run mode that allows the prosthesis wearer to limit the resistance in the joint to a minimum by switching over during activities such as cycling. The push button with bicycle symbol on the side of the joint now makes it even easier for the prosthesis wearer to switch to free run mode independently. The opposite push button allows the user to switch back to walking mode

The swing phase control is regulated by a powerful hydraulic damper. Thanks to the new design of the hydraulics, the damper can be conveniently adjusted from the front while the prosthesis wearer is standing.







The KINEGEN.stream is available in three versions. The following connections are possible proximally:



Adjustment core (3A2500)



M36 for 4-Arm Socket Adapter 3A2500/4



M36 for 3-Arm Socket Adapter (3A2500/3)











	3A860 KINEGEN.smart Brake Knee Joint with Lock	3A1000 KINEGEN.air Brake Knee Joint with Pneu- matic Swing Phase Control	3A900 KINEGEN.wave Polycentric Knee Joint with Lock - waterproof
max. body weight	125 kg	125 kg	125 kg
Mobility class (Walky):	1 to 2	1.5 to 3	1 to 2
Proximal Connection	Adjustment Core	Adjustment Core	Adjustment Core
Distal Connection	Tube Clamp Ø 30 mm	Tube Clamp Ø 34 mm	Tube Clamp Ø 30 mm
Knee Flexion Angle	145 °	145°	145°
Total Height	115 mm	249 mm	165 mm
Proximal Height	N/A	53 mm	N/A
Distal Height	N/A	196 mm	N/A
Overall Height of Joint Head	N/A	N/A	N/A
Effective Height	65 mm	185 mm	108 mm
Effective Proximal Height	27 mm	39 mm	25 mm
Effective Distal Height	38 mm	146 mm	83 mm
Weight	550 g	820 g	715 g













3A900/1 KINEGEN.wave Polycentric Knee Joint with Lock - waterproof

3A1800 KINEGEN air-active Polycentric Knee Joint with Pneumatic Swing Phase Control

3A2500 KINEGEN stream Polycentric Knee Joint with Hydraulic Swing Phase Control

3A2500/4 KINEGEN stream 3A2500/3 KINEGEN stream Polycentric Knee Joint with Hydraulic Swing Phase Control

Polycentric Knee Joint with Hydraulic Swing Phase Control

125 kg 125 kg 150 kg 150 kg 125 kg 1 to 2 2 to 3 3 to 4 3 to 4 3 to 4 M36 Thread Adjustment Core Adjustment Core M36 for M36 for 4-Arm Socket Adapter 3-Arm Socket Adapter Tube Clamp Ø 30 mm Tube Clamp Ø 30 mm Adjustment Core Adjustment Core Adjustment Core 145° 150° 136° 136° 136° 157 mm 228 mm 224 mm 217 mm 222 mm aforementioned 53 mm N/A N/A N/A aforementioned 196 mm N/A N/A N/A aforementioned N/A 36 mm 36 mm 30 mm 99 mm 185 mm 196 mm 209 mm 203 mm 16 mm 14 mm 23 mm 36 mm 30 mm 173 mm 173 mm 173 mm 83 mm 164 mm 700 g approx. 960 g 1110 g 1150 g 1120 g

## Accessories

For an all-round perfect fitting, all knee joints of the KINEGEN series have corresponding accessories.

The overview shows a selection of the most important articles.

You can find the complete selection and further information at www.streifeneder.com/op.



<b>3A1000 KINEGEN.air</b> Brake Knee Joint with Pneun Swing Phase Control	3A900 KINEGEN.wave natic Polycentric Knee Joint with Lock - waterproof	3A1800 KINEGEN.air-active Polycentric Knee Joint with Pneumatic Swing Phase Control and Integrated Extension Assist	3A2500 KINEGEN.stream Polycentric Knee Joint with Hydraulic Swing phase control
			×
×	×	×	
×		×	×
×		×	×





Streifeneder ortho.production GmbH Moosfeldstraße 10 82275 Emmering Germany Phone +49 8141 6106-0 Fax +49 8141 6106-70 export@streifeneder.de www.streifeneder.com/op



300W77/E 2021-09

This brochure is protected by copyright. The photomechanical reproduction of the entire brochure, even in extracts, is only permitted with the express permission of Streifeneder ortho, production GmbH. The stated dimensions and weights are non-binding guide values. During the delivery period, the manufacturer reserves the right to make changes in construction or form, deviations in color shade as well as changes in the scope of delivery, provided that the object of purchase is not significantly changed and the changes are reasonable for the buyer.

Changes may be made to the product after the editorial deadline for the brochure, 08.09.2021. The general terms and conditions apply, which you can view at www.streifeneder.com/op.

