ManuMax

96B1/...

Handgelenkbandage/Wrist support

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Instructions for use

Mode d'emploi

Istruzioni d'uso

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Brugsanvisning

Bruksanvisning

Käyttöohje

Bruksanvisning

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Upute za uporabu


Návod na používanie

Navodila za uporabo


Made in Germany


1

2

3

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EN

Purpose

ManuMax provides anatomical guidance for the affected wrist, thus relieving it during movement. The adjustable plastic rod and the textile strap provide mechanical stabilisation. Medical compressive active supports have a proprioceptive effect. They improve self-perception and can compensate for muscle imbalances. They reduce pain, support the removal of fluids from the tissue, thanks to their rheological effect, and provide relief during movement.

Here is how to put on your ManuMax correctly:

1. Open the velcro strap and slide your hand into the wrist support until the stabilising splint lies on the

inner surface of your hand. Put your thumb through the opening provided.

- 2.** Position the wrist support in such a way that the ring-shaped depression in the pressure pad encloses the wrist bone without pressure.
- 3.** Now close the velcro strap with sufficient tightness so that you feel an additional stabilisation of your wrist.

The correct fit of the wrist support is assured when it encloses the wrist tightly but does not apply a too strong (painful) pressure to the wrist.

The wrist support can be readjusted at any time by repeating steps 2 and 3 above. Rough, unfiled nails or sharp-edged rings etc. may damage the very fine knitted fabric.

Please note:

Only wear your wrist support when prescribed by your physician. Medical aids should only be dispensed by suitably trained staff of a medical products supplier. The effectiveness of the wrist support can only be realised if it is applied correctly. The wrist support is largely resistant to oils, ointments, skin moisture and other environmental influences. Certain combinations may negatively impact on the durability of the fabric. We therefore recommend that it gets examined regularly by the medical products supplier. Should the product become damaged, please consult your specialist supplier. Do not repair the wrist support yourself – this could compromise its quality and medical efficacy. For hygienic reasons, the wrist support is not suitable for treating more than one patient. The sewn-in textile label is important for the identification and traceability of the product. We therefore ask you not to remove the label under any circumstances.

Washing and care instructions

Please note the care instructions on the textile label. Please take the stabilising splint out of the wrist support and remove the strap. Do up the fasteners before washing! Wash the supports separately the first time (colour may run). Your support should be washed every day on a gentle or delicate wash cycle (30° C). Tip: a laundry net provides additional protection for the knitted fabric. Please do not use fabric softener! When hand washing, rinse well and do not wring. We recommend the gentle Juzo special detergent. You can reduce drying time by placing the support on a thick towel after washing, tightly rolling the towel up and firmly pressing out excess moisture. Then hang the support up loosely. Do not leave the product in the towel and do not dry it on a radiator or in the sun. The support is not suitable for tumble drying. Juzo supports are not suitable for dry cleaning. The removed stabilising splint can be cleaned with a damp cloth. Then slide it back into the pocket provided for it when the support is completely dry (with the bent spoon in front – the open side of the hollow section faces away from the hand surface). Finally velcro the strap onto the wrist support and it is ready for use again.

Material used

Details are provided on the label sewn into your wrist support.

Notes for the orthopaedic technician

The volar stabilising split is thermoplastically deformable at 180 °C. When shaping it, use suitable heat protection gloves to avoid injury.

Storage information and shelf life

Store in a dry place and protect from sunlight.

Supports have a maximum shelf life of 42 months. The usage period of the medical product is printed on the box with an hourglass symbol. The medical retailer or doctor prescribing the support will advise on exactly how long it can be used for. Body measurements may change based on the patient's condition and individual circumstances. We therefore recommend that medical specialists check body measurements regularly.

Indications

Distortions (sprains/strains), tendosynovitis (tendovaginitis), arthrosis (degeneration of the wrist), arthritis (inflammation of the wrist), post-traumatic irritation of the wrist

Contraindications

If you have any of the following conditions, you should only wear the support after consulting with your doctor: Lymph drainage disorders and non-specific soft tissue swelling outside the area of application, sensory disturbances and impaired circulation in parts of the body to be treated, skin conditions, skin irritation and injuries in the area to be treated; open wounds in the treatment area must be covered with sterile material. Liability for failure to comply with the contraindications cannot be accepted by the Julius Zorn GmbH Group.

Side effects

There are no known side effects when used as directed. However, if negative alterations (such as for example skin irritation) should occur while the use of our products has been prescribed, please go to your doctor or your specialist dealer immediately. If an incompatibility should be known against one or several ingredients of this product, please check back with your doctor before use. If your symptoms worsen during the time you are wearing the product, please take off the support and consult a medical specialist immediately. The manufacturer is not liable for damages/injuries which have been caused by improper handling or misuse.

In the event of complaints related to this product, such as damage to the material or problems with the fit, please contact the medical retailer directly. Only serious incidents that may lead to a major deterioration of the patient's medical condition or to death should be reported to the manufacturer and the relevant authority in the Member State. Serious incidents are defined

in Article 2 (65) Regulation (EU) 2017/745 (MDR).

Disposal

As a rule, the product is to be destroyed via the conventional disposal route after the end of use. Please observe the applicable local/national regulations.

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Vertrieb/Distributor: Streifeneder · ortho.production GmbH
Hersteller/Manufacturer: Julius Zorn GmbH
Berro AG · Austrasse 7 · 9055 Bühler · Suisse
UK responsible person: Juzo UK Ltd