


Sturdy and flexible at the same time: Aramid Fibre Fabric 95P29

The tearproof aramid fibre fabric is perfectly suitable to manufacture partially flexible areas in forefoot- or thigh orthoses. Due to its high durability (E-module), it is very flexible compared to carbon fibres and therefore ideal for layer construction to arm laminate compounds. The material can be optimally cut with the electric cutter 80C 10.

Characteristics

- high durability, impact- and abrasion resistance
- high weight-related tensile strength
- good shock absorption performance
- very high flexibility
- store away from UV-rays





Streifeneder
ortho.production GmbH
Moosfeldstrasse 10
82275 Emmering
Germany
T +49 8141 6106-0
F +49 8141 6106-70
export@streifeneder.de
www.streifeneder.com/op

300W144/E 2017-06