

Active daily living with the Komet
The new therapeutic shoe in outdoor-design



Comfortable wound care



Modern, flexible, high quality – with its new outdoor-design, the latest forefoot-relief shoe Komet is the perfect everyday companion for patients. The subtle and restrained look is provided by the especially low rocker sole, which makes a balance shoe unnecessary for the therapy. Despite its low-profile concept, the sole is secure and stable due to the slip-proof profile.

Post-surgical forefoot relief shoe

With its wide, square sole shape in the forefoot area, the Komet is optimally suitable for foot care after surgery techniques with Kirschner-wire. Further indications are peripheral arterial occlusive diseases, wounds in the forefoot area due to diabetes mellitus as well as independent wounds in the forefoot area.

Protection and hygiene

Due to its high-quality materials, the Komet offers maximal wearing comfort. Reliable protection from dirt, water and cold temperatures is provided together with the forefoot protection cap 217R12, which can be ordered separately. A special feature is the anti-bacterially equipped cover fabric of the padded sole: an ammonium-silane compound limits the growth of certain bacterial strains and thus supports hygienic wound care.



Sensory loss due to diabetic polyneuropathy (PNP), i.e. occurring sensations such as pressure, temperature, pain of the feet cannot be noticed in time. We recommend to regularly check feet and shoes.

All features at one glance



Shoe size	Inner shoe length	Inner shoe width	Item-no.
35 – 36	approx. 242 mm	approx. 103 mm	217R1/36
37 – 38	approx. 255 mm	approx. 106 mm	217R1/38
39 – 40	approx. 268 mm	approx. 109 mm	217R1/40
41 – 42	approx. 282 mm	approx. 112 mm	217R1/42
43 – 44	approx. 295 mm	approx. 115 mm	217R1/44
45 – 46	approx. 308 mm	approx. 118 mm	217R1/46
47 – 48	approx. 322 mm	approx. 121 mm	217R1/48
49 – 50	approx. 335 mm	approx. 124 mm	217R1/50

Streifeneder
ortho.production GmbH
Moosfeldstrasse 10
82275 Emmering
Germany
T +49 8141 6106-0
F +49 8141 6106-70
export@streifeneder.de
www.streifeneder.com/op



300W146/E 2018-03

This brochure is protected by copyright.

Reproduction of the complete brochure, even in abstracts, in any form or by any means, is only permitted with explicit written authorisation of the Streifeneder ortho.production GmbH.

The indicated measurements and weights are guide values without obligation. During the delivery period, changes in design or form, deviations in colour as well as modifications of the score of delivery shall remain reserved on part of the manufacturer, as long as the purchased product is not substantially altered and as long as the changes are reasonable for the buyer.

Products are subject to change without prior notice after press date of this brochure, 2018-03-28.

Our general terms and conditions apply; they can be viewed anytime at www.streifeneder.com/op.