


Fast and efficient for plaster models: Drying Hose 99P23

As a perfect aid in silicone technology, the drying hose offers the possibility of drying plaster models almost completely in just 3 – 4 minutes with compressed air. This requires the use of porous moulding plaster 5E1.

Characteristics

- dries porous plaster in no time
- hose end must be sealed in the model
- stabilisation of the hose by attachment to the plaster rod
- drying hose can be reused
- can be cut to the required length





Streifeneder
ortho.production GmbH
Moosfeldstrasse 10
82275 Emmering
Germany



T +49 8141 6106-0
F +49 8141 6106-70
export@streifeneder.de
www.streifeneder.com/op

300W154/E

