


Comfort Pads for orthotic closures 145R14

More wearing comfort for users: the pre-finished pads are pulled through the loops into the orthotic closure. They can be placed individually, protect optimally against pressure spots and can be flexibly trimmed.

Characteristics:

- available in two sizes and colour (blue, black)
- material: skin side polyamide (PA); outside polyester (PES) und polyvinylchlorid (PVC) (blue), polyester (PES) and polyamide (PA) (black)
- suitable for trimming at the embossing
- for strap widths: 25 and 30 mm (M); 38 and 50 mm (L)
- disinfectable with common alcoholic disinfectants
- PU = 10 pcs.





Streifeneder
ortho.production GmbH
Moosfeldstrasse 10
82275 Emmering
Germany



T +49 8141 6106-0
F +49 8141 6106-70
export@streifeneder.de
www.streifeneder.com/op

300W157/E