


High quality Magnetic Closure for orthoses 44P20/25

Simple operation, high quality workmanship: the magnetic closure Fidlock[®], made of synthetic material and stainless steel is easy to close and open with just one hand. The system is suitable as spare part for PeroSupport.tec 197R4 as well as for the construction of individual orthoses.

Characteristics

- Weight: 22 g; weight load capacity: 70 kg
- Dimensions: 74 x 33 x 13 mm
- Strap passage: 2 mm
- for strap width 25 mm
- PU: 1 pc.





Streifeneder
ortho.production GmbH
Moosfeldstraße 10
82275 Emmering
Germany



T +49 8141 6106-0
F +49 8141 6106-70
export@streifeneder.de
www.streifeneder.com/op

300W160/E