

Reinforces support with every step

Vacuum adapter 7P50

The new active vacuum adapter 7P50
for a comfortable and secure connection from
the residual limb to the prosthesis:

It intensifies the hold through negative pressure inside
the socket. Lifting movements are thereby reduced.



Designed for a strong connection

The 7P50 adapter actively creates vacuum while walking:

After only 8-10 steps, negative pressure of up to -0.5 atm is achieved and maintained by active pumping.



7P50 At a glance

- made of aluminium
- Pipe diameter 30 mm
- up to 100 kg body weight
- proximal connection adjustment core support
- including accessories

Description	Item. no:
Vacuum adapter	7P50

