

The new

KINEGEN.comfort

comfortable and versatile

With its modern, aesthetic **design** combined with **sophisticated technology**, the KINEGEN.comfort fits seamlessly into the KINEGEN series.

It combines comfort with safety:
due to **three-phase pneumatic swing phase control, load-dependent brake** and **rear-mounted articulated axle**.

The KINEGEN.comfort is suitable for a **wide range of activities** and medium walking speeds. It enables a harmonious gait pattern and safe walking – even on ramps.



Advantages for prosthetists and patients

- convenient brake and pneumatic **adjustment from the front**
- **integrated mechanics** are easy on the cosmetic foam
- **harmonious gait pattern** thanks to three-phase pneumatic swing phase control
- **wide range of activities**
- **high stance phase stability**
- **safe walking and standing**, even on ramps



Technical data

designation	effective installation height	weight	max. body weight	activity level	knee flexion angle	item-no.
KINEGEN.comfort	180 mm	847 g	125 kg	1,5-3	145°	3A1010